

SPECIALTY **SPRINGS**



GOODMAN TORQUING SPRING

Gain full control of your treatment with the Goodman Torquing Spring™. This spring can be placed to gain lingual or labial root torque on any tooth. Root alignment is important in preventing relapse and nothing moves roots faster, with better accuracy and control, than the Goodman Torquing Spring. Available in two sizes to fit any tooth.



TO LEARN MORE:

Use a code reader on your smart phone to scan this code



GOODMAN TORQUING SPRING 5 PACK

SIZE	ITEM NUMBER
Small	60.40.506.00001
Large	60.40.507.00001

Large springs are for maxillary centrals, cuspids, and possible maxillary laterals.

Small springs are for mandibular incisors and small maxillary laterals.

Springs fit .016 x .022, .017 x .025, .016 x .025, .018 x .025, and .019 x .025.

ANTERIOR TORQUING SPRING

60.80.800.05064

5 PACK

The Anterior Torquing Spring generates more positive torque in the upper incisors. The device ends with its hooks distally to each upper central bracket, and the center circle (coil) pressed into the archwire.

